


### Art-Line™: Arterial Catheter Sets



Code	Description	Box quantity
HP-1822	18Ga X 8cm Arterial Catheter	10
HP-1841	18Ga X 12cm Arterial Catheter	10
HP-2000	20Ga X 4.5cm Arterial Catheter	10
HP-2022	20Ga X 8cm Arterial Catheter	10

#### Intent of use

Arterial catheter is generally use to provide continuous arterial pressure measurements. Biometrix arterial catheter sets present all in one solution for performing arterial catheterization.

#### Technical specification

- HP-18XX:
  - 18Ga Radiopaque arterial catheter
  - 0.035" x 50cm marked guidewire, J flexible tip



- 18Ga x 7cm introducer needle
- 5ml syringe
- Scalpel
- HP-20XX
  - 20Ga Radiopaque arterial catheter
  - 0.021" x 20cm straight guidewire
  - 20Ga x 4cm introducer needle
  - 3ml syringe

### **Quality and regulatory information**

The quality system is in compliance with ISO 13485:2012, ISO 9001:2008 and the Medical Devices Directive 93/42/EEC ANNEX II.

\*The information given in the present leaflet is a subject for change. Please contact the company representative for a confirmation.